eyepower *red*

Deep red light to promote eye health

INSTRUCTION GUIDE

English

Deutsche

Francais

Espagnol

WELCOME

Welcome to eyepower and thank you for purchasing our unique deep red light eyewear.

eyepower *red* is based on research carried out by University College London (UCL), Department of Ophthalmology and other universities around the world and is designed to help promote healthy eyes as you age in a safe and convenient way.


SAFETY

Deep red light at 670nm occurs naturally in sunlight and has been approved as safe by the Ethics Committees of University College London and the world-renowned Moorfields Eye Hospital.

The power and wavelength of eyepower *red* have been independently tested and verified that they conform to the same values as those used by UCL in their research studies.

DESIGN

To make eyepower *red* as affordable as possible we have created a 'one-size-fits-all' design. We recommend that the device is used when sitting or lying down to give you the most effective and comfortable experience.

On the underside of the device is a charging cable socket and an on button which also acts as an indicator light.

HYGIENE

To clean your eyepower *red* simply wipe it down with an anti-bacterial wipe. The device is NOT waterproof and should not be exposed to liquids.

HOW TO USE THE DEVICE

For best results please observe the following points:

Duration

Studies have shown that a 3 minute session of deep red light is sufficient to have a positive biological effect on the cells of your retina which will help keep your eyes healthy.

To start the session simply press the On button on the base. The red light will then come on at low brightness for 3 seconds before going to full brightness for 3 minutes. At the end of the 3 minutes the light will turn off automatically and the indicator light on the base will flash green 5 times to show a session has been completed successfully.

Frequency

The beneficial effect of a single 3 minute session will last for around 7 days. Therefore, as a minimum it is recommended to have at least one session every week.

Additional sessions may provide some additional benefits, particularly for specific eye conditions. However, the overall outcome will vary between individuals and will also depend on the type and severity of the eye condition. Please seek professional medical advice if you have any concerns.

If you have normal healthy eyes then 2 or 3 sessions a week may well provide additional benefit but there is currently insufficient data to quantify the precise advantage of more frequent use. Deep red light has very low energy and so more frequent sessions are not a risk to healthy eyes.

Time Of Day

Based on your body's natural circadian rhythm, your eyes are at their most receptive to red light in the morning. Using red light in the afternoon or evening has little or no beneficial effect.

Eyes Open Or Closed

If you find the light of *eyepower red* too bright, it can be used with your eyes closed. Your eyelids are almost transparent to light at 670nm so the same benefit will be achieved.

Colour Changes

After the session has started it is most likely that you will start to see the red light change to different colours. It is quite common to see green, orange, white, yellow or even a mix of all these colours. This is completely normal and is due to the way the eye reacts to the red light.

After Image

At the end of the session you will experience an after image for a short period of time. This is most often seen as green but might also show objects with a strong red tinge. This is perfectly normal and will fade after a few minutes as your eyes adjust back to normal lighting levels.

CHARGING THE BATTERY

Your *eyepower red* uses an internal rechargeable battery to power the device. As the electronics are always 'on' there is a continual drain on the battery, so depending on usage, you will need to recharge it every 2 to 3 weeks. However, for convenience you may want to give the battery a weekly boost rather than waiting for it to run low. The indicator light will flash blue depending on the battery and charging state as outlined below:

Connecting The Charging Cable Correctly

A red USB 'C' type cable is supplied with the *eyepower red*. The larger end is plugged in to a USB charger plug, while the smaller 'C' plug is pushed firmly in to the socket on the base of the *eyepower red*. The indicator light will flash blue depending on the battery and charging state. Please ensure that you push the 'C' plug firmly into the base - it is quite a tight fit and will not move around once it is fully inserted.

eyepower red showing the charging cable fully inserted.

The cable is held firmly in place and should not loose when touched.

Not Fully Charged After Being Connected For A Few Hours

If your *eyepower red* is not showing a constant blue light (fully charged) after a few hours on charge then please check the following:

- 1) Make sure the cable is connected correctly as outlined above.
- 2) If the plug is connected correctly and no constant blue light after 5 or 6 hours, then your USB charger plug may not deliver enough current to charge the device fully.

If this is the case, then try swapping to a different charger plug. Most modern mobile phone charger plugs, such as Apple or Samsung, work well.

Indicator Lights

The indicator light on the base of the *eyepower red* will light up blue as shown below to indicate various charging states:

If the indicator light flashes red, green and blue then it indicates that the battery is out of charge. To clear this, plug your *eyepower red* in to charge and give one short press on the ON button. After a few seconds the light will give the normal 3 blue flashes to show it is charging.

You may find it more convenient not to wait until the low charge light flashes and just charge your *eyepower red* once a week or every two weeks anyway, it will not do it any harm.

SESSION REMINDER

eyepower red has a visual reminder function using the indicator light that flashes on a weekly basis to remind you that another red light session is due. The Session Reminder is turned off by default. If you wish to use the session reminder, please visit the [eyepower web site support page](http://www.eyepower.co.uk/support) at www.eyepower.co.uk/support for instructions on how to use it.

TECHNICAL DETAILS

Light Source:	Dual Array LED SMD
Wavelength:	670 - 675nm Near Infrared
Light Output:	8.0 - 8.5mW/cm ²
LED Lifetime:	> 20,000 Hours
Flicker Frequency:	No flicker
Battery:	120mA/Hr Single Cell LiPo Battery
Charging:	USB cable with 'C' Type Connector
Warranty:	2 Years From Shipment Date

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HELP & SUPPORT

If you have a query relating to the use of **eyepower *red*** with a pre-existing eye condition please refer to your clinical professional or optician.

eyepower *red* is not intended as a medical device and our staff are not qualified to provide a medical opinion.

For more information on using your **eyepower *red***, including information on the Session Reminder, please visit:

www.eyepower.co.uk/support

For a large print copy of this instruction guide please download from:

www.eyepower.co.uk/guide

If you cannot find an answer to your query in the support section then please email details to:

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Light Power Health Ltd